

German **ACCURACY**



Rapid STAT

AUSTRALIAN STANDARD

MEETS AS 4760 TARGET CUT-OFFS

COMPONENT	LEVEL
AMPHETAMINE (AMP)	50NG/ML
COCAINE (COC)	50NG/ML
METHAMPHETAMINE (MET)	50NG/ML
OPIATES (OPI)	50NG/ML
CANNABIS (THC PARENT)	25NG/ML

CMM Technology™ can ensure supply of stock from our temperature controlled facility, so you are not caught short.

Call today: **1300 79 70 30**

SCREENING FOR ORAL FLUID

Rapid STAT is a point of collection drug test for the simultaneous detection of five drugs of abuse in oral fluid specimen within 8-12 minutes.

It is a lateral flow immunoassay in which each drug of abuse analyte is represented by a separate line in the test window of the device. The test system consists of an untreated sample collector, a buffer bottle and a test device.

Rapid STAT is intended for use in the qualitative detection of the following drugs of abuse in human oral fluid. Amphetamine, Methamphetamine, Cocaine, Opiates, THC.



 **CMM**
Technology™

